

The highest demands are placed on packaging in the medical technology and healthcare sector today, from maximum protection against mechanical damage and effective protection against contamination to double sterile barrier packaging for maximum security against the penetration of germs.

Whenever you need high-quality, practical and safe protective packaging – you have found the right place. More information online.



PROTECTIVE COVER



TUBES



BOXES



CASES



STORAGE SYSTEMS



TipCover

TipCover is tiny – but big in protection

03/2026

rose plastic
medical packaging GmbH

Rupolzer Straße 30
D-88138 Hergensweiler/Lindau
Tel. +49 8388 92345-0
info@rose-medipack.de



Should we make a big thing out of our tiny TipCover?

Like so often in life, small things can make a big difference. Our TipCover is such a small thing. A tiny one, actually. Its 8 sizes range from 28 to 36 mm length.

Yet TipCover is perfectly equipped to do what it is intended to do. It protects the sensitive tips of your medical devices such as implants or surgical instruments from damage. We won't make a big thing out of this. It's a tiny thing, after all. Although one with big benefits.



TipCover –
protects sensitive
tips from people.
And protects sensitive
people from tips.



PRODUCED IN CLEANROOM

Produced in cleanroom, it's suitable for sterile as well as non-sterile applications.

Evenly important, it protects product users from injuries that sharp edges or tool tips can cause so easily. And when customer can't get hurt, that won't hurt your customer relations either.

100 % SILICONE-FREE

Tip Cover is silicone-free. No silicone can get in contact with patients.



Perfect tip
protection

TipCover

Perfect tip protection for sensitive tools or implants. 8 sizes cover a wide range of requirements and dimensions.



The proof of the suitability of the packaging in combination with the medical device to be packed must be provided by the manufacturer.

THE ADVANTAGES AT A GLANCE:

- Ideal for the protection of tips and edges of surgical instruments
- for sterile and non-sterile applications
- secure loading and removal
- easy handling
- safe transportation
- 8 sizes available for various diameter ranges

Product Range

Article	Outer dimensions W x H x D	Suitable for	
TCV0	0.5×28/1	10×28×5mm	∅ 0,5 – 2,5mm
TCV1	2.5×28/1	9,3×28,3×4,8mm	∅ 2,5 – 4,0mm
TCV2	4×28/1	12,5×28,4×7,5mm	∅ 4,0 – 7,0mm
TCV3	7×28/1	17,0×28,5×12mm	∅ 7,0 – 10,5mm
TCV4	10.5×28/1	23,3×28,7×16mm	∅ 10,5 – 14,5mm
TCV5	14.5×36/1	29,6×35,8×22mm	∅ 14,5 – 18,0mm
TCV6	19×2×28/1	22,4×28×12mm	2,5×8,0 – 3,5×16,0mm
TCV7	25×3×28/1	28,4×28×12mm	2,5×13,0 – 4,0×22,0mm